


Instructions for Use








---

RE-FLEX SHOCK<sup>™</sup>

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
DA	Brugsanvisning	25
SV	Bruksanvisning	29
NL	Gebruiksaanwijzing	33
PT	Instruções de Utilização	37
JA	取扱説明書	41
ZH	中文说明书	44

## **A COMPLETE RE-FLEX SHOCK FOOT SYSTEM COMPRISES THE FOLLOWING:**

- Re-Flex foot module
- Maintenance free shock module
- A replaceable external front spring attached to the shock module
- For category 9 a replaceable titanium coil spring located inside the shock module
- Cosmetic foot cover

## **OPTIONS:**

Re-Flex foot modules are available either as Standard (**Figure 1**) or Tall (**Figure 2**).

**WARNING:** Avoid exposure to fire, excessive heat or cold. This may damage the shock or change its function.

## **PYRAMID ADAPTER**

Use Loctite 243 and torque to 15Nm (11ft-lbs).

## **FLEX-FOOT SOCK**

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area.

## **SHOE HORN**

When removing and installing the foot cover use a Flex-Foot shoehorn to avoid damage to foot and/or foot cover.

## **ALIGNMENT RECOMMENDATIONS**

To ensure optimum function of the foot and pylon, use the following guidelines:

### **1. Bench Alignment**

To introduce the appropriate socket flexion and heel height: Divide the foot length into 3 equal portions (**Figure 3**). The load line should fall at the junction of the posterior and middle third portions of the foot (**Figure 3**).

### **2. Adjusting the Shock Module**

To adjust the shock module loosen lower pin on the spring (**Figure 4**) and the three hex screws (**Figure 8**) on the attachment clamp. Rotate the shock housing until correct height is achieved (**Figure 5**). One full rotation results in 8mm change in height, height can be set in half-turn steps of 4mm.

The spring (shown in grey colour in figure 6) can be rotated 180° in relation to the threaded tube (shown in white in figure 6) so the spring can always be positioned at the front when adjusting the height.

**WARNING:** Ensure there is full engagement of the clamp and shock module. Less than full engagement will void the warranty. Leave a small gap between the collar and the clamp.

### 3. **Securing the Shock Module**

To secure the shock module fasten the lower pin on the spring (**Figure 4**) applying 10Nm with a torque wrench and tighten the 3 hex screws (**Figure 8**). Start with the top hex screw and work down by applying 3Nm with a torque wrench. Perform this sequence twice, since the upper hex screw will loosen in the first run. Place the bellow in the groove on the threaded tube and then tighten the collar over the bellow. (**Figure 7**)

**WARNING:** Over tightening the hex screws can strip threads and cause shock module to bind.

**WARNING:** Make sure the collar covers the edge of the bellow all around.

### 4. **Applying a Permanent Shock Module Setting**

To permanently fix the shock module, remove the hex screws one by one; apply Loctite 243. When finished, re-apply torque of 3Nm on all hex screws. Perform this sequence twice, since the upper hex screw will loosen in the first run. Failing to apply Loctite 243 to the hex screws may result in the shock coming loose.

Also remove the pin and apply Loctite 243 on threads and tighten with 10Nm torque.

## **DYNAMIC ALIGNMENT**

### **Heel Wedges**

Heel wedges influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The small, medium and large wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics. The wedges can be trimmed using sharp scissors to customize stiffness.

### **Temporary Wedge Placement (Figure 10)**

- Place the wedge in the angle of heel and foot module.
- Secure in position with tape wrapped around the foot module.

### **Permanent Wedge Placement**

- Roughen the upper and lower surface of the heel with abrasive paper.
- Apply adhesive on the lower side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module. (**Figure 10**)

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

### **Heel Divider**

The heel divider should be placed in the last third of the split heel part. To secure its position a drop of instant adhesive can be used on one side.

## CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	326-366
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	N/A
High Impact Level	2	3	4	5	6	7	8	9	N/A	N/A

## CLEANING

The product needs to be cleaned after contact with water, sand or dirt. The filters need to be cleaned with air pressure. (Figure 9)

## MAINTENANCE

No component maintenance.

The filters need to be cleaned with air pressure (Figure 9).

**Change spring:** Remove the top and bottom pins using 4mm “hex-key” and change the spring (Figure 4). Apply Loctite 243 on threads and tighten pins with 10Nm torque.

**Change bushing in spring:** Remove spring (see above) and press out the bushings carefully. Press new bushings in and re-attach the spring (see above).

## LIABILITY


The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guideline.

## COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been adapted adequately on the

basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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