

Instructions for Use

RE-FLEX ROTATE™


		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	20
DA	Brugsanvisning	24
SV	Bruksanvisning	28
NL	Gebruiksaanwijzing	32
PT	Instruções de Utilização	36
PL	Instrukcja użytkowania	40
CS	Návod k použití	44
RU	Инструкция по использованию	48
JA	取扱説明書	52
ZH	中文说明书	56


1


2


3


4

5

6

ENGLISH

A COMPLETE RE-FLEX ROTATE FOOT SYSTEM COMPRISES THE FOLLOWING:

- Re-Flex foot module
- Maintenance free shock module with rotational compliance
- A cosmetic foot cover (ordered separately).

OPTIONS:

Re-Flex Rotate is available with either a male or female pyramid connector.

Re-Flex foot modules are available either as Standard (**Figure 1**) or Tall (**Figure 2**).

WARNING: Avoid exposure to fire, excessive heat or cold. This may damage the shock or change its function.

PYRAMID ADAPTER

Use Loctite 243 and torque to 15Nm (11ft-lbs).

FLEX-FOOT SOCK

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area.

SHOE HORN

When removing and installing the foot cover use a Flex-Foot shoehorn to avoid damage to foot and/or foot cover.

ALIGNMENT RECOMMENDATIONS

To ensure optimum function of the foot and pylon, use the following guidelines:

1. Bench Alignment

To introduce the appropriate socket flexion and heel height:

- Divide the foot length into 3 equal portions (**Figure 4**).
- The load line should fall at the junction of the posterior and middle third portions of the foot (**Figure 4**).

2. Adjusting the Shock Module (**Figure 1**)

To adjust the shock module, loosen the three hex screws on the attachment clamp (**Figure 3**) and rotate the foot until the correct height and/or rotation is achieved. One full rotation results in 8mm change in height. For finer adjustments, rotate in 90° turns:

Warning: Ensure there is full engagement of the clamp and shock module. Less than full engagement will void the warranty.

Turn:	Build height change
-90°	2mm lift
-180°	4mm lift
-270°	6mm lift
-360°	8mm lift
+90°	2mm lowering
+180°	4mm lowering
+270°	6mm lowering
+360°	8mm lowering

3. Securing the Shock Module

To secure the shock module tighten the 3 hex screws (**Figure 3b**). Start with the top hex screw and work down by applying 3Nm with a torque wrench. Perform this sequence twice, since the upper hex screws will loosen in the first run.

Warning: Overtightening the hex screws can strip threads and cause shock module to bind.

4. Applying a Permanent Shock Module Setting

To permanently fix the shock module, remove the hex screws one by one; apply Loctite 243 and 3Nm with a torque wrench. When finished, re-apply torque of 3Nm on all hex screws. Failing to apply Loctite 243 to the hex screws may result in the shock coming loose.

DYNAMIC ALIGNMENT

Heel Wedges

Heel wedges influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The small, medium and large wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics. The wedges can be trimmed using sharp scissors to customize stiffness.

Temporary Wedge Placement (Figure 5)

- Place the wedge in the angle of heel and foot module.
- Secure in position with tape wrapped around the foot module

Permanent Wedge Placement

- Roughen the upper and lower surface of the heel with abrasive paper.
- Apply adhesive on the lower side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module (**Figure 6**).

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

Heel Divider

The heel divider should be placed in the last third of the split heel part. To secure its position a drop of instant adhesive can be used on one side.

PROTECTIVE FILM ON PYRAMID

The protective film protects the pyramid from scratching. Remove the film when the fitting is complete.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Low Impact Level	1	1	2	3	4	5	6	7	8
Moderate Impact Level	1	2	3	4	5	6	7	8	N/A
High Impact Level	2	3	4	5	6	7	8	N/A	N/A

MAINTAINANCE & CLEANING

Note: The device should always be used with spectra socks and foot covers to minimize dust and dirt from entering between carbon blades. It is important to fix the spectra sock around pyramid location in order to prevent noise issues coming from sand entering between the plates. The device is weather-proof but not resistant to corrosion. Therefore, the device should not come into contact with salt water or clorinated water. The device is not designed to be used under extreme conditions like diving or jumping into water. If for some reason the foot components get wet parts should be dried with a lint free cloth. Should the device come into contact with salt water, chlorinated water, dust or dirt it must be rinsed with fresh water immediately.

Note: The stiffness of the shock module can increase in temperatures below -15 °C (+5 °F).

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal

body mass. In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur					
Category	Weight (Kg)	Lable text			
1	52	ISO 10328 -	P3	52	kg
2	59	ISO 10328 -	P3	59	kg
3	68	ISO 10328 -	P3	68	kg
4	77	ISO 10328 -	P4	77	kg
5	88	ISO 10328 -	P4	88	kg
6	100	ISO 10328 -	P5	100	kg
7	116	ISO 10328 -	P6	116	kg
8	130	ISO 10328 -	P7	130	kg

- This manual is intended for use by a certified prosthetist.

WARRANTY

Össur provides a warranty of 36 months for this foot module and 2 months for the foot cover

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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