









Instructions for Use

VARI-FLEX[®] MODULAR

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	25
DA	Brugsanvisning	30
SV	Bruksanvisning	35
NL	Gebruiksaanwijzing	40
PT	Instruções de Utilização	45
JA	取扱説明書	50

ENGLISH

VARI-FLEX® MODULAR FOOT SYSTEM (Figure 1)

- Minimal Clearance 182mm (7")
- Maximum user weight 166 kg (365 lbs)
- Available in sizes 22-30 cm

BENCH ALIGNMENT (Figure 2+3)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram (the mark on the brim of the foot cover matches the one third portion).

CUTTING FOOT MODULE (Figure 4)

Before cutting check for enough clearance

- Measure for length and add 25mm (1") for T Connector, 38mm (1 1/2") when using a pyramid.
- **Warning:** The T Connector/pyramid should only be placed where it makes full contact with the carbon blade.
- Clamp foot module in vise. Use a sharp fine tooth hacksaw blade.
- Do not cut below the cut limit sticker line the lower cutline for Attachment.
- Sand away any sharp edges.
- Secure foot to pyramid or T using C-clamp
- Check height of prosthesis on patient while standing between parallel bars.
- Cut to final length.

DRILLING FOOT MODULE ATTACHMENT HOLES

- Back-up foot module when drilling using wood block or discarded carbon fiber section and clamp securely.
- Use drill guide template to locate holes (Figure 5).
- Use sharp drill bits. Bullet nose or solid carbide tipped drill bits suggested.
- Remove template and drill to appropriate hole size.

Drill size:		
Cat 1-6	8,5mm	21/64"
Cat 7-9	10mm	25/64"

- Drill speed 1200-1800 RPM.
- Feed drill slowly with light intermittent pressure - do not force.
- Foot is now ready for assembly.

USING MODULAR ALIGNMENT ADAPTER

- Attach the cutting adapter to the foot module (Figure 9).
- Make sure grit is between adapter and foot module.
- Connect the foot and the cutting adapter to the prosthesis.

- Adjust the cutting adapter to the correct height by sliding it on the blade.
- Clamp the cutting adapter to the foot module by tightening the screws to 22Nm.
- Disconnect the foot module, with the cutting adapter clamped on to it, from the prosthesis.
- Drill appropriate holes according to "Drilling Foot Module Attachment Holes" instructions above through the drill bushings instead of the guide template. Note that the drill bushings can only accommodate 8.5mm drills.
- Therefore a second drilling with a 10.5 drill is required for the larger categories.
- Mark or measure correct length of foot module.
- Detach cutting adapter from the foot module.
- Cut foot module according to "Cutting Foot Module (**Figure 4**)" instructions above.
- Attach pyramid to foot module.
- Connect the foot module with the pyramid attached to the user.

DYNAMIC ALIGNMENT

The heel stores energy at heel strike, then slowly releasing it to assist tibial progression to midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	Special Order
High Impact Level	2	3	4	5	6	7	8	9	Special Order	Special Order

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

Actions

- Add polyurethane heel wedge.
- Shift socket anterior (or foot posterior).

Heel too hard

Symptoms

- Rapid heel to toe movement.
- At initial contact the amputee has poor control of his prosthesis.

- Minimal energy return feeling.
- Knee may become unstable.

Actions

- Replace heel wedge with a smaller one.
- Shift socket posterior (or foot anterior).

Heel Wedges

The wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The small, medium and large wedges can be interchanged to customize the stiffness

and achieve the desired functional characteristics. The wedges can be trimmed using sharp scissors to customize stiffness.

Temporary Wedge Placement

- Cut the wedge to the width of the foot module.
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module with the curved surface of wedge towards the foot module.
- Secure in position with tape wrapped around the foot module (Figure 6).

Permanent Wedge Placement (Figure 7)

- Apply adhesive on the lower side of the wedge only (1).
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module (Figure 8).

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

Heel Assembly

If it is necessary to replace the heel the heel component use Loctite 410 and torque to following specifications:

Category	Bolt Size	Heel Bolt Torque
1-6	1/4"	14 Nm (10ft-lbs)
7-9	5/16"	28 Nm (21ft-lbs)

Heel Divider

The heel divider should be placed in the last third of the splitted heel part. To secure its position a drop of instant adhesive can be used on one side.

FLEX-FOOT SOCK

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area and secured proximally with the cable tie provided.

SHOE HORN (Figure 10)

When removing and installing the cover use shoehorn to avoid damage to foot or cover.

ENVIRONMENTAL CONDITIONS

- Operating temperature: -15 °C (+5 °F) to 50 °C (122 °F).
- Operating humidity: 10-95% relative humidity
- Shipping and storage temperature: -20 °C (-4 °F) to 70 °C (158 °F).

Note: The device should always be used with spectra socks and foot covers to minimize dust and dirt from entering between carbon blades. It is important to fix the spectra sock around pyramid location in order to prevent noise issues coming from sand entering between the plates.

WATERPROOF IF USED WITH CARBON TUBE/DIRECTLY LAMINATED ONTO THE SOCKET

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

WEATHERPROOF IF USED WITH PYRAMID/T CONNECTOR


A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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