






Instructions for Use

---

ELATION<sup>®</sup>

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	14
DA	Brugsanvisning	16
SV	Bruksanvisning	18
EL	Οδηγίες χρήσης	20
NL	Gebruiksaanwijzing	23
PT	Instruções de Utilização	25
JA	取扱説明書	28
ZH	中文说明书	30
KO	사용 설명서	32

# ENGLISH

---

## **ELATION® FOOT SYSTEM (Figure 1)**

Build Height:

- Male Pyramid 87mm (3 7/16")
- Adjustable heel height: 0-5cm (0 - 2")
- Weight Limit: 100kg (220lbs)
- 4 interchangeable heel bumpers
- 1 double durometer forefoot bumper
- For low to moderate impact levels
- Sizes: 22 – 28cm

## **BENCH ALIGNMENT (Figure 2 + 3)**

- Adjust to appropriate heel height (using the shoe). Do not exceed the recommended heel height.
- Introduce appropriate socket angles; flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.

## **DYNAMIC ALIGNMENT**

I Coronal Plane

*Check for:*

- Relative positioning of socket and foot.
- Contact with ground should be even.

II Sagittal Plane

*Check for:*

- Smoothness of gait.
- Socket to foot position.
- Heel stiffness.

## **PUSH BUTTON MECHANISM**

- Press the black button to unlock the ankle.\*
- Adjust the ankle to the appropriate heel height by using body weight on the toe or heel.
- Pylon or shin section should be perpendicular to the floor when the foot is in the shoe.
- Sit down when pressing the white button.
- Lift the foot off the ground to reach the button.
- Press the white button to lock the ankle.

\*The foot cannot be used in that position for weight bearing

**NOTE:** Using an extension to the push button can result in an accidental release.

## **HEEL AND FOREFOOT BUMPERS (Figure 4)**

**One forefoot bumper, 4 heel bumpers available**

To exchange bumpers in order to increase or decrease the stiffness of the heel:

## ATTENTION!

Use a sharp hex key. Special care needs to be taken when untightening/tightening bolts for sizes 26-28 as the hex slot is shallow. If the hex slot is worn out replace with a new bolt (**Figure 5**).

- Remove mounting bolts with hex key.
- Compress the heel bumper in order to expose the forefoot bumper (a).
- Remove bumper.
- Position foot in maximum dorsiflexion to expose and remove heel bumper (b).
- Replace with desired bumper ensuring the bumper is fully seated.
- Reassemble in reverse order.
- Apply Loctite 243 to mounting bolts, reassemble ankle to foot and torque M6 bolts to 15Nm on sizes 22-25 and M8 bolts to 21Nm on sizes 26-28.

Heel Bumper Selection Guide		
Weight	Durometer	Heel
50 kg - less (110 lbs - less)	Soft	Green
51 kg - 65 kg (111 lbs - 143 lbs)	Medium	Black
65 kg - 80 kg (143 lbs - 175 lbs)	Firm	Red
80 kg - 100 kg (175 lbs - 220 lbs)	Extra Firm	Orange

**NOTE:** A Dual Density Shin Fairing FCX0002 is available. It allows the range of motion in the ankle section to be maximized without compromising the integrity of the foam where it is attached to the foot cover. The soft section also facilitates access to the push button on Elation.

## LIABILITY


The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guideline.

## COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks.

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
**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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