













Instructions for Use

BALANCE™ KNEE

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	24
NO	Bruksanvisning	29
DA	Brugsanvisning	33
SV	Bruksanvisning	37
EL	Οδηγίες χρήσης	41
FI	Käyttöohjeet	46
NL	Gebruiksaanwijzing	50
PT	Instruções de Utilização	55
PL	Instrukcja użytkowania	60
CZ	Návod k použití	65
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	73
JA	取扱説明書	78
ZH	中文说明书	82
KO	사용 설명서	86

NUMBERING IN FIGURES

1. Alignment reference line
2. Socket
3. Pivot Axis
4. Front links
5. Chassis
6. Housing
7. Back link
8. Lock
9. Lock Disable Cap
10. Release Cord
11. Handle assembly
12. Plastic tube
13. Metal fasteners
14. Crimp
15. Screw
16. Screw Insert
17. Collar washer

FUNCTIONS AND BENEFITS

- Adult polycentric knee with geometric stability
- Excellent stability at heel strike and mid-stance yet easy to release with toe load
- Increased toe clearance in swing phase
- Easy adjustment of stability setting, extension assist and swing friction
- Available with a manual lock function
- Up to 180° of knee flexion(165° with the manual lock)
- 24 months limited warranty

INDICATIONS FOR USE

- The Balance Knee is intended exclusively for lower limb prosthetic fittings
- Adult amputee up to 125kg (275 lbs)
- Single speed limited community ambulator

CLEARANCE (Figure 1)

Specific instructions regarding Manual Locking Feature

Alignment notice

The Balance Knee joint is delivered already adjusted for ideal locking position. Care must be taken when adjusting knee stability not to compromise the function of the locking mechanism.

ALIGNMENT

Bench alignment

- Assemble prosthetic foot, pylon and knee.
- Establish knee centre height, measured from the pivot axis.
- If pylon is cut too short, do not compensate by putting a spacer into the knee tube adapter. Doing so may compromise strength and voids the warranty (**Figure 2**). Use a new Össur pylon cut to correct length for maximum user safety.
- Tighten tube clamp bolt to 15 Nm (11 ft-lbs)(**Figure 3**) with a 5mm hex key. Do not overtighten.

- Place appropriate heel height under heel or align with user's shoe.
- Attach socket with selected adapter combination to knee.
- Ensure flexion/extension and abduction/adduction angles are as determined in user assessment.
- Make sure the adapter does NOT touch the back link or front links during full flexion (**Figure 4**). This can cause damage to the knee and will void the warranty. Contact between socket and chassis or pylon during full flexion is acceptable (**Figure 5**).
- The alignment reference line, from bisection of the socket on ischial level, should pass 0-10 mm in front of the pivot axis of knee and through posterior 1/3 mark on the Flex-Foot cover (**Figure 6**).

Lanyard mechanism (Figure 7)

- The release cord should run along the lateral or anterior side of the prosthesis.
- Laminate screw insert and plastic tube into the socket at a desired location. For test sockets, a hole can be drilled through the socket wall followed by pre-heating the thread insert and pressing it in place from inside the socket.
- Attach Balance Knee to the socket.
- Fasten the handle assembly to socket with screw and collar washer in place.
- Thread release cord through plastic tube(if used) and metal fasteners on each end.
- Thread release cord through unlocking split ring/handle and back through metal fastener.
- Shorten the release cord to the length that suits the user and make sure that the release cord is slightly tensioned.
- Crimp fasteners by pliers or pinch the center with a punch and mallet.
- Ensure that the release cord is not impaired

Static alignment

- For safety, please make initial adjustment with user standing between parallel bars!
- Fit prosthesis and check for correct length.
- Make sure foot is flat on the floor, knee is in neutral position and socket angles are correct.
- When user shifts weight onto the prosthesis, the locking function should be activated or the knee should be in stable position. Modify alignment if necessary.

PRODUCT SPECIFIC TRAINING AND ADJUSTMENTS

Instructing new users is essential to achieve a successful rehabilitation.

DYNAMIC ALIGNMENT

For safety, please make initial adjustments with the user between parallel bars!

Explain function of knee stability. It is important the user understands how the knee stability functions before walking.

- When the locking mechanism is activated the flexion of the Balance Knee Joint is blocked.
- When the locking mechanism is deactivated the Balance Knee Joint is stable from heel strike to midstance.
 - Knee flexion is initiated during roll-over of the foot.

- Ask the user to actively release the knee from stable position by moving the weight over the toe simulating roll-over and then back to foot flat.
- Roll over can also be simulated by pelvic rotation and it is important that the user is aware of this, to ensure stability when turning.

Locking Mechanism settings

- Use the lock disable cap to temporarily activate or deactivate the locking mechanism(**Figure 8**). To deactivate the lock put the lock disable cap on the lock and to activate the lock remove the lock disable cap.
- To deactivate the locking mechanism permanently: remove the knee cap (**Figure 9**). Remove the two split rings (**Figure 10**). Use 2.0 mm hex key to remove the lock pivot pin (**Figure 10**). Remove both lock and springs (**Figure 11**). Put the knee cap back on.

LEVEL GROUND WALKING

Adjust the stance release setting if needed:

Stance release setting (Figure 12):

- Use 4 mm hex key to adjust the stance release position, this setting must be adjusted without weight bearing to prevent damage.
- Clockwise: Earlier stance release and initiation of flexion becomes easier.
- Counter-clockwise: Later stance release and increased stability.

When a suitable knee release setting has been found, the extension assist can be adjusted to the user's self selected walking speed.

Extension assist setting (Figure 13):

- Insert a 4 mm hex key in the round holes of the grey spring housing from the rear to adjust spring tension.
- Right turn: Increases extension assist.
- Left turn: Decreases extension assist.

When preferred extension assist setting has been achieved the constant friction can be adjusted to eliminate terminal impact if necessary. **Please note that excessive constant friction may cause insufficient extension speed.**

Constant friction setting (Figure 14):

- Use 4 mm hex key to adjust the constant friction, be careful not to overtighten, max torque is 3 Nm.
- Clockwise: Increases friction.
- Counter-clockwise: Decreases friction.

As confidence increases, have the user walk outside of the parallel bars. Make sure the user can change directions and turn around without using pelvic rotation on the prosthetic side, as this may cause inadvertent knee joint release.

SITTING DOWN

- Place feet evenly (or prosthetic foot slightly anterior) and shift weight onto forefoot using hip flexors. Knee will initiate flexion.
- Use sound limb to accept weight for sitting

- Users finding it difficult to sit down can also benefit from using pelvic rotation to initiate flexion of the knee instead of hip flexors.

STAIR DESCENT

For safety reasons it is recommended that Balance Knee users descent stairs with a fully extended knee.

FOLLOW UP CHECK

- Reassess adjustments after a couple of weeks
- It may be necessary to change adjustments after user has gained experience with the knee.

CARE AND MAINTENANCE

- The Balance Knee is designed for a maintenance free service life.
- Excessive exposure to water/moisture or submerging the knee may however significantly reduce the service life and will void the warranty.
- If cleaning is needed, wipe the knee with soft cloth moistened with small amount of general purpose oil (3 in 1) or sewing machine oil.
- DO NOT use solvent to clean the knee. Bearings and seals may be damaged.
- DO NOT use compressed air to clean knee. Air forces pollutants into bearings and may cause malfunctions and wear.

CAUTION! Avoid placing hands or fingers near moving joints.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to 3 million load cycles. Depending on the amputee's activity this corresponds to a duration of use of 3 to 5 years. We recommend carrying out regular yearly safety checks

ISO 10328 - P6 - 125 kg *)

*) Body mass limit not to be exceed!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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
EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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