





Instructions for Use





---

TOTAL KNEE<sup>®</sup> 2000  
Geometric Locking System

		3
EN	INSTRUCTIONS FOR USE	6
DE	GEBRAUCHSANWEISUNG	13
FR	NOTICE D'UTILISATION	22
ES	INSTRUCCIONES PARA EL USO	30
IT	ISTRUZIONI PER L'USO	39
DA	BRUGSANVISNING	47
SV	BRUKSANVISNING	54
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	61
NL	GEBRUIKSAANWIJZING	70
PT	INSTRUÇÕES DE UTILIZAÇÃO	79
PL	INSTRUKCJA UŻYTKOWANIA	87
CS	NÁVOD K POUŽITÍ	96
TR	KULLANIM TALIMATLARI	104
RU	ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ	112
JA	取扱説明書	121
ZH	中文说明书	128
KO	사용 설명서	135


11

2052


2053


2047

12


2075R


2075L

13

27

14

28

28

# ENGLISH

---

## NUMBERING IN FIGURES

1. Alignment reference line
2. Socket
3. Knee centre reference
4. Front link
5. Balance unit
6. Pivot axis (alignment axis)
7. Hydraulic housing
8. Back link
9. Bottom link
10. Chassis
11. Bumper
12. Shim
13. Bumper cover
14. Extension promoter
15. Valve F (Flexion)
16. Valve H (Heel Rise / Initial Flexion)
17. Valve E (Extension)
18. Spring cup
19. Spring
20. Ferrule
21. Cable
22. Retaining Ring
23. Spring cup seat
24. Knee cap – small
25. Knee cap – large
26. Shin ferrule
27. Bearings (lubrication points)
28. Release points

## FUNCTIONS AND BENEFITS

- Increased knee stability and cadence variation.
- Increased toe clearance in swing phase.
- Lightweight design – slim profile for enhanced cosmetic finish.
- 160° knee flexion – suitable for sitting, kneeling, cycling and other activities where greater knee flexion is desired.
- 24-month warranty.

## INDICATIONS FOR USE

- Adult amputee up to 100kg (220lbs) patient weight.
- The Total Knee is exclusively used for lower limb prosthetic fittings.
- Suitable for an unlimited community ambulator who has potential to walk with varying cadence.

## CLEARANCE (Figure 1)

## ALIGNMENT AND ADJUSTMENT (Figure 2)

## ADAPTER OPTIONS (Figure 3)

- Prong adapter (A-114040 or A-114030) for knee disarticulation or long transfemoral residual limbs.
- Female pyramid (2026) for medium to short transfemoral residual limbs.

- Euro 4-Hole (2055) for medium/short transfemoral residual limbs, in combination with 4-Hole Socket Adapter (A-122100 or A-125100).
- Male pyramid (2025 or 2025T) for medium/short transfemoral residual limbs or hip disarticulation.

Important: Always use original Total Knee/Össur adapters. Fully engage threads. Align pinch bolt 90° anteriorly for male pyramid, prong adapters and Euro 4-hole adapter, but 45° for female pyramid (Figure 3). Threads are cut for maximum strength in this position.

Adapters from other manufacturer may not fully engage all threads. They may compromise strength and void the warranty. Never tighten adapters against the hydraulic head.

When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the knee.

Important: Contact between socket/adapter and proximal 2/3 of back link must be avoided when fully flexing knee (Figure 4). Contact in this area can damage knee and void the warranty. Contact between socket and back link below balance unit is acceptable.

If pylon is too short for changed setup, do not compensate by putting a spacer into the knee tube adapter as this may compromise strength and void the warranty (Figure 5). Use a new Össur pylon cut to correct length for maximum patient safety.

### **BENCH ALIGNMENT (Figure 2)**

- Assemble prosthetic foot, pylon and knee.
- Establish knee centre height, measured from the knee center reference(3). Tighten tube clamp to 15 Nm **(Figure 6)**.
- Place appropriate heel height under heel or fit patient's shoe.
- Establish vertical pylon for proper knee function. Alignment reference line (1) from pivot axis (6) should pass through posterior 1/3 mark on Flex-Foot® cover (or see other foot system instructions).
- Attach socket with selected adapter combination to knee. Tighten to torque as specified in adapter IFU. Ensure flexion/extension and abduction/adduction angles are as determined in patient assessment. Alignment reference line (1) from bisection of socket on ischial level (2) should pass through pivot axis (6) of knee and posterior 1/3 mark on Flex-Foot® cover.
- Make sure total length of prosthesis is correct.

NOTE: If alignment reference line (1) is too far anterior to pivot axis(6), more effort may be required to initiate knee flexion.

### **BUMPER SELECTION (Figure 7)**

The bumper determines the amount of stance flexion.

- Carefully insert flat head screwdriver or similar under bottom of bumper cover and remove cover (Figure 8).
- Insert small flathead screwdriver under bumper and lift upwards.
- Remove bumper (leave the shims under it in place).
- Install appropriate bumper according to patient weight.

<b>Patient Weight</b>	<b>Colour</b>	<b>Type</b>
35-50 kg (80-110 lbs)	Green	Soft

<i><b>Patient Weight</b></i>	<i><b>Colour</b></i>	<i><b>Type</b></i>
50-70 kg (110-155 lbs)	Black	Regular
70-90 kg (155-200 lbs)	Red	Firm
90-100 kg (200-220 lbs)	Yellow	Extra Firm

- Re-install the bumper cover by hooking the groove under the chassis front and pushing the two pins into the chassis holes.
- Active patients may require a stiffer bumper
- Patients who have previously worn a prosthesis without a stance-flexion feature may initially prefer a stiffer bumper.
- We suggest trying a softer bumper after two weeks of trial. This helps the patient develop a security with the knee. Doing this will also help make the transition to the stance-flexion characteristic more acceptable.

### **STATIC ALIGNMENT**

For safety, please make initial adjustments with patient standing between parallel bars!

- Fit prosthesis and check for correct length.
  - Make sure foot is flat on the floor, knee is in neutral position and socket angles are correct.
  - When the patient shifts weight onto the prosthesis, the geometric lock is activated and the bumper should be slightly compressed. Modify alignment if necessary.
  - Explain function of stance flexion/geometric lock to patient. Ask patient to activate it by stepping forward with prosthesis and transferring weight onto heel.
  - Patient should become accustomed to activating geometric lock and resulting stance flexion movement.
  - Knee cannot collapse as long as geometric lock is activated.
  - Preset swing flexion resistance, Valve F (**15, Figure 9**)
- a. Sit patient on chair. Extend prosthetic knee. Let it fall into flexion.
  - b. Increase resistance of Valve F until a small bump observed at 60° flexion angle.

### **DYNAMIC ALIGNMENT**

The model 2000 knee provides smooth and fluid swing control during normal ambulation and allows changes in cadence.

For safety, please make initial adjustments with patient standing between parallel bars!

#### **Valve adjustment – F, H, E (Figure 9)**

<i><b>Valve</b></i>	<i><b>Factory setting</b></i>	<i><b>Function</b></i>
Valve F	3/4 open	Affects swing flexion resistance from 60°-160° (fast walking only)
Valve H	3/4 open	Affects swing flexion resistance from 0°-60° (slower walking)
Valve E	1 1/2 open	Affects whole swing extension resistance, 160°-0°

Use small flathead screwdriver for adjustment.

- If patient is walking fast and you observe excessive heel rise, increase resistance of Valve F until heel rise has normalised.
- Use Valve H to smooth flexion resistance changes between initial flexion and 60°.
- Increase swing extension resistance with Valve E until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension as Geometric Lock is activated prior to heel contact. Excessive swing extension resistance will prevent knee from moving into full extension and eliminate activation of Geometric Lock.
- Fine tune Valves F and E until a smooth and secure gait pattern is achieved for slow and fast walking speeds.

### **HEEL RISE / INITIAL FLEXION – VALVE H (Figure 9 & 10)**

Generally it is not required to adjust Valve H, unless the patient's walking speed is slow and excessive heel rise is observed.

Valve H can also be used to fine tune the resistance between initial flexion and 60°.

To adjust Valve H remove the „sticker“ (Figure 10).

**CAUTION!** Over tightening valves or flexing knee with all valves completely closed may damage valves.

### **Bumper adjustment**

- Try different bumpers and select the most comfortable one.
- Final bumper selection may vary from bumper indicated for patient's weight according to table.
- Soft bumpers increase the locking moment strength of the Geometric Lock and increase stance flexion motion.
- Hard bumpers decrease the locking moment of the Geometric Lock and decrease stance flexion motion.
- Foot selection influences bumper selection.
  - a. For soft heels, choose a softer bumper.
  - b. For hard heels, choose a hard bumper.

### **Shim adjustment (Figure 7)**

Adding shims changes the release point of the Geometric Lock allowing for swing phase knee flexion. This affects knee stability and ease of walking and sitting. The shims are thin plastic plates which are installed underneath the bumper to adjust its height. A change in height affects the knee's release point.

- More shims: easier release into swing flexion and decreased stance flexion
- Less or no shims: increased safety, increased stance flexion and delayed release into swing flexion

**NOTE:** There are two different types of shim. The thicker ones are 0.76mm thick and clear, the thinner ones are 0.38mm thick and grey. We do not recommend using more shims than those two combined. Adding more shims can eliminate the Geometric Lock and affect stability of the knee during stance.

Gradually increase shim thickness until the user feels premature release of the knee. Reduce shims until the user finds a balance between security and release point. Please make initial adjustments inside parallel bars.

**NOTE:** If no shims are used, the knee will exhibit some flexion/extension play during stance. If the patient finds this motion unacceptable, this can be reduced by adding shims.

### **EXTENSION PROMOTER (Figure 7)**

The extension promoter (14) reduces excessive heel rise and promotes faster knee extension. **It is set at minimum when shipped from Össur.**

- Flex the knee to access extension promoter adjustment screw.
- To reduce excessive heel rise, especially when walking fast, turn the extension promoter adjustment screw clockwise with a large flathead screwdriver.

### **EXTENSION ASSIST (Figure 10)**

The extension assist is an **OPTIONAL** feature that helps get the knee back to full extension. We recommend use of the extension assist only when the amputee is unable to extend the knee to engage the geometric lock.

**NOTE:** Prior to installing the extension assist, try reducing Valve E resistance to help achieve full extension.

- Beneficial for users having difficulty achieving full extension.
- Can help during fast walking.
- Brings knee to full extension consistently.
- Places positive locking bias on Geometric Lock to enhance stability.

Available with three different springs (yellow=regular, blue=firm, red = extra firm).

### **INSTALLATION (Figure 10)**

- Remove bumper cover, bumper and shims (Figure 7 & 8)
- Install spring cup (18) into hole inside 30mm tube clamp of knee (23)
- Firmly seat spring cup by gently tapping with rubber mallet
- Place ferrule (20) over cable (21)
- Place spring (19) over cable (21) and ferrule (20) (start with the yellow spring)
- Insert cable (21) and spring (19) into spring cup (18), ball end of cable first
- Slide ball end of cable through hole in knee (23) and insert ball into the slot on bottom link (9) using small screwdriver or needle nose pliers.
- Reinstall the selected shims (12), bumper( 11) and bumper cover (13).
- Let patient walk, and adjust the swing extension with Valve E. If sufficient extension is not obtained with a fully open Valve E, try a stiffer spring.
- Readjust swing flexion resistance after installation of extension assist (Valves F and H)

### **FOLLOW-UP CHECK**

- Reassess adjustments after a couple of weeks
- It may be necessary to change shim and bumper adjustments after user has gained some experience with knee
- Increased experience, security and voluntary control may require change to harder bumper and/or increased shim thickness.

### **COSMETIC OPTIONS (Figure 11 & 12)**

We recommend a discontinuous cosmetic solution which does not interfere with knee function. Continuous cosmetic solutions are also possible.

## Discontinuous Cosmetics

Ref. No	Description	Notes
2047	Shin ferrule	Protects unit from foam cover - maintains space for movement within foam cover
2052	Knee cap - small	Use to prevent tearing at the knee
2053	Knee cap - large	Use for patients who kneel often
2075R	Foam cover - regular	Accommodates shin ferrule
2075L	Foam cover - large	Accommodates shin ferrule

## MAINTENANCE

The Total Knee 2000 is designed for low-maintenance. It is warranted against defective materials and workmanship for 24-months. We recommend:

- Total Knee 2000 be inspected every six months for signs of unusual wear, by a prosthetist familiar with this product.
- Stance flexion bumper must be replaced if cracked or worn, depending on patient's activity level.
- If knee is subjected to excessive moisture or corrosive environment, it's recommended to clean and lubricate the knee frequently.

## CLEANING

- Wipe knee with soft cloth moistened with small amount of kerosene. DO NOT use solvent stronger than kerosene. If kerosene is not available, use general purpose oil (3in1) or sewing machine oil.
- DO NOT dip knee or pour solvent over knee. Bearings and seals may be damaged.
- DO NOT use compressed air to clean knee. Air forces pollutants into bearings and may cause malfunctions and wear.

## LUBRICATION (Figure 13)

Three sets of roller bearings are visible. Place a few drops of machine oil on the roller bearings (general purpose oil (3in1) or sewing machine oil can also be used). Move the knee several times and wipe off excessive oil with a soft cloth. No other parts need external lubrication.

## PRODUCT-SPECIFIC TRAINING

Instructing new users is essential to achieve a successful rehabilitation. For safety, initial adjustments and gait training should be done with patient standing between parallel bars.

## LEVEL GROUND WALKING

- Explain to the user the function of stance flexion/geometric lock. Patients who previously wore a knee without stance flexion function may fear collapse of knee during loading response. They must be trained to allow the knee to move into stance flexion and not to prevent it by strong hip extension
- Explain the stance release mechanism to the patient. Patient's centre of mass must be anterior to the pivot axis of knee. In this position, very little hip flexor effort is required to initiate knee flexion. Knee flexion will be initiated through rotation of the hip as weight is transferred onto the ball of the foot during terminal stance.

- Allow patient to stand in the parallel bars with their weight over the ball of the foot. Initiate knee flexion with pelvic rotation and slight hip flexion. Repeat several times.
- Start walking within parallel bars. As confidence increases, start walking outside of the parallel bars.

## SITTING DOWN

- Have the patient move their centre of mass forward over the toe
- Place feet evenly (or prosthetic foot slightly posterior), shift weight anterior onto forefoot. Knee will release
- Use sound limb to accept weight for sitting

## STAIR DESCENT

Walking step over step down stairs with Total Knee requires much practice and should only be attempted by trained Total Knee users. Begin practicing on lowest step and use handrail.

- Place entire foot on step.
- Move centre of mass over forefoot as for sitting. Knee flexion will be initiated. Knee flexion will be rapid. Use caution!
- Increase number of steps as patient becomes confident.

Not all users will be able to learn step over step stair descent.

## MANUAL FLEXION

When sitting with fully extended Total Knee it may be necessary to flex knee manually. Press two points shown in **Figure 14** to release knee into flexion. Be careful not to pinch fingers in the linkages. (Patient may also reach out and pull or tap toe of prosthesis to release).

**CAUTION!** Avoid placing hands or fingers near moving joints.

## LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## COMPLIANCE

This component has been tested according to ISO 10328 standard to 3 million load cycles. Depending on the amputee's activity this corresponds to a duration of use of 3 to 5 years. We recommend carrying out regular yearly safety checks

**ISO 10328 - P5 - 100 kg \*)**

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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