







Instructions for Use

RHEO KNEE[®] XC

3

EN		Instructions for Use	4
DE		Gebrauchsanweisung	16
FR		Notice d'utilisation	28
ES		Instrucciones para el uso	40
IT		Istruzioni per l'uso	52
NO		Bruksanvisning	64
DA		Brugsanvisning	76
SV		Bruksanvisning	88
EL		Οδηγίες Χρήσης	100
FI		Käyttöohjeet	112
NL		Gebruiksaanwijzing	124
PT		Instruções de Utilização	136
PL		Instrukcja użytkowania	148
CZ		Návod k použití	160
TR		Kullanım Talimatları	172
RU		Инструкция по использованию	184
JA		取扱説明書	197
ZH		中文说明书	208
KO		사용 설명서	219

DESCRIPTION OF SYMBOLS

Attention – See instructions for use

On / Off

Type B Applied Part (IEC60601-1)

CE label with Notified Body identification number

Non-ionizing radiation

Manufacturer and date of manufacturing

Direct current

Consult Instructions for use

This marking on the product, packaging, accessories or literature indicates that the product contains electronic components and/or batteries that should not be disposed of in regular waste at the end of its usable life. To prevent possible harm to the environment or human health from uncontrolled waste disposal users are requested to separate these items from other types of waste and recycle them responsibly to support sustainable reuse of material resources. Users should contact their local government office for information on how these items can be recycled or disposed of in an environmentally sound manner.

To protect natural resources and to promote material reuse, please separate batteries and electrical components from other types of waste and recycle them through your local, free electronic parts return system.

Battery status

Serial number

Protected against solid foreign objects of 2.5 mm and greater and protected against splashing water.

The RHEO KNEE XC is referred to as the "*device*" in the following. This document provides you important information on the indication for use, alignment, setup and handling of the device. This manual is intended for use by a certified prosthetist and the user of the device. The device may only be fitted by qualified personnel authorized by Össur after completing the corresponding training.

PRODUCT DESCRIPTION (Figure 1)

The device is a microprocessor prosthetic knee joint and consists of the following components:

1. Pyramid adapter
2. Actuator
3. Extension lock
4. Power button
5. Led light (indicator for battery)
6. Charge plug
7. Rechargeable battery pack
8. Pyramid adapter

Make sure that the end user of the device knows how the device functions and how to operate with it. Please inform the end user about all safety precautions before he starts using the device!

FUNCTION

The device automatically adjusts its behavior to activities such as standing, walking, and descending stairs/ramps.

Standing

Whenever stance is detected the device locks and provides the user with full support. This allows the user to take small steps more naturally without having a locked knee feeling.

Walking

The device supports the user whenever ground contact is detected, releases as swing initiates and provides swing phase resistance allowing the user to achieve comfortable and symmetric transitions at a wide range of speeds. The device reduces the perceived exertion and thereby allows for walking longer distance.

Sitting down

As the user leans back to sit down the device responds to proportional loading to support sitting down with a natural movement.

Descending stairs and ramps

The device flexes and supports the user while lowering to the next step. It is recommended to lean slightly backwards as this facilitates the prosthetic knee bending and, in case of loss of balance, prevents the user from falling forward.

Ascending stairs

Stair ascent mode allows the user to ascent stairs step over step. The mode can be activated leading with the sound side or the prosthetic side. Placing the prosthetic foot onto the first step with the prosthetic knee in a slightly bent position and transferring weight onto the foot triggers stair ascent function. Alternatively, the user can place the sound side foot on the first step, properly flexing the hip on the prosthetic side to swing the device through, place it on the next step and transfer weight onto the foot. After the first step the device adapts to stair climbing and allows the user to ascent naturally.

CAUTION: Always use the handrail when operating the Rheo Knee in stairs. Particular caution is required when carrying children in stairs.

Cycling

The cycling function is automatically activated after completion of 2 crank rotations and will remove knee resistance allowing for efficient cycling. The device indicates cycling function with one long beep. When stepping off the bike, the user needs to fully extend the knee to return to normal knee function. Return to normal knee function is confirmed by one short beep.

Running

The device provides seamless transition from walking to running. Knee function automatically adjusts to running whenever running gait is recognized.

MEDICAL INDICATION

The device is to be used exclusively for the exoprosthetic fitting of amputations of the lower limb.

- Unilateral transfemoral amputation
- Bi lateral transfemoral amputation
- Unilateral knee disarticulation amputation
- Bi lateral knee disarticulation amputation

INDICATIONS FOR USE

- Activity level: Moderate to high
- The device is weatherproof, meaning it withstands exposure to humid and rainy weather conditions without damage or loss of function

Impact Level	Maximum weight limit
Moderate	136 kg (300 lbs)
High	110 kg (243 lbs)

CONTRAINDICATION FOR USE

Any conditions contrary outside or beyond to the indication for use and/or guidelines contained in this document.

SAFETY

WARNING: Statement that alerts the user to the possibility of serious injury, or serious adverse reactions associated with the use or misuse of the device.

CAUTION: Statement that alerts the user to the possibility of a problem with the device associated with its use or misuse.

PRECAUTIONS: Please follow the safety precautions as failure to do so can lead to malfunction of the device and risk injury to the user. Prosthetists should make sure their clients are aware of, and understand, the safety precautions.

GENERAL SAFETY INSTRUCTIONS

It is highly recommended that the user operates with the power off prior to setting swing and stance controls, but always use caution when walking with the power off.

1. For maximum safety use handrail when descending stairs and ramps. Particular caution is required when carrying children down stairs.
2. Improper handling or adjustment of the device may cause the stance control to malfunction, which can increase the risk of fall.
3. Make sure that the user is able to recognize the auditory and/or vibratory signal for power loss.
4. When malfunctions occur or the power reaches critical levels, auditory and/or vibratory warnings alert the user. If warning signals are detected, check if the knee can provide stance control. Engage the extension lock and proceed to walk with caution.
5. If the device loses power, the knee loses its ability to provide stance support. In such cases it is recommended to manually engage the integrated extension lock to lock the knee in full extension. The knee allows full weight bearing when locked.
6. Avoid any impact to the socket, frame, adapter contact, or the interface panel on the back of the device. Damage to this panel or any of its components may cause the product to malfunction (**Figure 4**).
7. The device may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

8. Do not attempt any repairs (no matter how minor) on the device. This may damage the knee, cause malfunction or failure, and will void the warranty.
9. Never operate this product if it is not working properly, if it has been dropped or damaged. In such cases return the product to your healthcare provider for examination and repair.
10. The device is weatherproof but is not designed to be exposed to large amount of water or other fluids, such as rinsing, showering or submersion.
If the device has been submitted to any of those conditions it must be sent in for service to Össur.
11. The device is not resistant to corrosion and hence the device should not come into contact with salt or chlorinated water. If the device has been submitted to any of those conditions it must be sent in for service to Össur.
12. The device is not resistant to ingress of solid particles smaller than 2.5 mm (1/10 in) and therefore the device should not come in contact with sand or be used in dusty environment. If there is any sign of impaired motion of the knee joint it must be sent in for service to Össur.

CHARGING THE DEVICE

WARNING: Make sure the user knows the risk of using damaged power supply, battery charger or plug can lead to electric shock.

Immediately replace damaged power supply units, plugs or battery charger.

CAUTION: Make sure that the user knows how to check the charge level before use and charge the device if required. The Össur Logic app can be used to track battery status.

Before fitting the device, charge the unit for 3 hours to ensure full charge.

HOW TO CHARGE THE DEVICE

- Insert the charger cable into the charging port on the user interface panel (**Figure 1**).
- It will take approximately 1.5 hours to charge the device battery to 90%, and 3 hours to fully charge it (100%).
- The device should be charged daily for optimal performance.
- It is not possible to overcharge the device.
- The device can be charged regardless of whether it is turned ON or OFF

BATTERY WARNINGS (Figure 3)

When battery is low a repeating warning pulse will appear: A short sound repeats every 15 seconds and orange light flashes on the interface panel. At this point the device should be recharged. How much time is left on the battery depends on activity (approximately 30 minutes of walking on level ground or 5 min of stair descent).

When battery is critically low, a buzzing, vibrating and a flashing red light will signal an imminent shut down (after 5 seconds).

BATTERY CAPACITY

Rechargeable Lithium Ion Battery (1950 mAh/28.9 Wh): Fully charged battery is sufficient for 48–72 hours of continuous use, depending on activity.

CAUTION:

- Follow local ordinances or regulations for proper disposal of device, accessories and packaging. The battery should be disposed of in accordance with local laws and regulations. It should not be incinerated.
- Because the device is battery operated, caution should be taken about excessive mechanical shocks or vibrations.

TURNING THE DEVICE ON AND OFF

HOW TO TURN THE DEVICE ON (Figure 1)

To turn the device on press the ON/OFF button. You will hear a 2 beep confirmation.

When the device is on, you will see a pulsing green light on the light indicator beside the ON/OFF button on the user interface panel.

HOW TO TURN THE DEVICE OFF (Figure 1)

To turn the device off press the ON/OFF button and hold it down until you hear 4 short beeps and 1 long beep. When the device is off, the light indicator on the user interface panel will also turn off.

EXTENSION LOCK

See (Figure 6).

CALIBRATING THE DEVICE

REQUIREMENTS FOR CALIBRATION

In order to calibrate the device, the accompanying software is required (RHEOLOGIC software for PCs, Össur Logic app for iOS devices)

The accompanying software technical manuals give instructions on how to connect to the software and modify adjustable parameters for the device.

HOW TO CALIBRATE THE DEVICE

1. Before starting make sure there is NO LOAD on the knee and DO NOT connect the socket to the knee.
2. Make sure the knee is FULLY extended.
3. Connect to the knee using the accompanying software.
4. Open the “Calibrate” tab in the accompanying software and press the “Calibrate” button.
5. When calibration is complete the device will provide a sound / vibration signal.

The following values are acceptable calibration ranges:

- Toe load (-2.5 to 2.5)
- Heel load (-2.5 to 2.5)

WARNING: Improper calibration (for example if load is applied on the knee or socket is connected) affects the knee's ability to perform optimally and may compromise the user's safety.

RECOMMENDED ALIGNMENT OF THE DEVICE

RECOMMENDED COMPONENT SELECTION

Please refer to the list of recommended components for the device in the Össur catalog.

BENCH ALIGNMENT

1. Connect the device to the foot (with the foot cover on) using corresponding Össur adapters / connectors.
2. Establish knee center height.
3. Mount the socket to the device using corresponding Össur adapters / connectors.
4. The alignment reference line should pass through a sagittal plane bisection of the socket, the knee center and the posterior 1/3 of the foot. (Figure 2)

NOTE: If stability modifications are desired after patient assessment; the reference line can be shifted up to 5 mm (3/16") anterior or 5 mm (3/16") posterior to the knee center.

CAUTION: As the alignment reference line is positioned more posterior to knee center, the user will experience more stance flexion during loading response. Additional voluntary control will be required to maintain knee stability when the device is powered OFF.

PRECAUTION: At maximum flexion, maintain minimum distance of 3 mm (1/8") between the device and the socket. At maximum flexion where contact cannot be avoided due to socket volume, make sure that the contact is on the frame, directly below the user interface panel. For maximum pressure distribution maintain a flat and cushioned contact surface (Figure 2).

STATIC ALIGNMENT

1. Ask the user to stand with equal weight on each foot and make sure that the pelvis is level (Figure 5).

2. Evaluate sagittal and coronal alignment for correct position of the alignment reference line.
3. The alignment reference line should pass through a sagittal plane bisection of the socket, the knee center and the posterior 1/3 of the foot module.

NOTE: If stability modifications are desired after patient assessment; the reference line can be shifted up to 5 mm (3/16") anterior or 5 mm (3/16") posterior to the knee center.

DYNAMIC ALIGNMENT

The best way to acquire a dynamic alignment is by having a symmetric gait pattern.

The initial dynamic assessment is important for getting the user familiar with relying on the prosthesis. Össur recommends that this assessment is completed at the parallel bars. Complete general dynamic alignment at the same time.

1. Evaluate the user's level of voluntary control by having the user walk with the device powered off. The user should have the ability to maintain knee stability when the power is off.
2. Use hip extension to maintain stability into stance.
3. Have the user experience the response, the roll-over, of the foot.
4. Load the heel and allow the knee to flex, thus experience the resistance provided by the knee.
 - a. Control the rate of the knee flexion with hip extension.
 - b. Perform rapid loading of the prosthesis with use of the parallel bars.
5. Train the user to maintain normal step length on the sound side.
6. Ask the user to sit down using the stance resistance of the device.

CAUTION: Errors that can occur during alignment and adjustment of the device, can result in malfunction of the device. This can cause the user to fall.

COMPLETING SETUP

SOCKET/ FRAME CONTACT POINTS

Care should be taken when configuring the socket / frame and adapters so that contact only occurs at the posterior apex of the device.

DO NOT remove the devices protective cover. This will immediately void the warranty.

ACCOMPANYING SOFTWARE

Prosthetists have the option of using RHEOLOGIC Workbench for Windows PC and Össur Logic in expert mode for Apple iOS devices to set up and adjust parameters in the knee.

The end user can use Össur Logic in the default user mode to monitor basic knee parameters such as battery status and step count.

For more information please refer to the instructions for use for the individual software products.

SPECIFICATIONS

STRUCTURAL SPECIFICATION

- Net weight: 1.6 kg (3.5 lbs)
- Frame construction: Aluminum
- Knee flexion available: 120 degrees
- Clearance and dimensions (**Figure 2**)

SOFTWARE COMPATIBILITY

RHEOLOGIC™ software for PCs or Össur Logic software for iOS devices.

PRODUCT SPECIFICATION INFORMATION

The device is tested and certified to comply with the IEC/EN60601-1, standard of electrical safety of medical devices and IEC/EN60601-1-2, electromagnetic compatibility for medical electrical devices and ISO 10328. The company fulfills the requirements of ISO13485, MDD 93/42/EEC, and the device is CE marked accordingly.

ENVIRONMENTAL CONDITIONS

- Operating temperature: -10°C to 40°C (14°F to 104°F)
- Operating humidity: 0% - 100% relative humidity
- Storage temperature: 0°C to 45°C (32°F to 113°F)
- Storage humidity: 10% - 90% relative humidity (non-condensing)

CAUTION: If the unit is either used or stored in cold environment (below 0°C) for a period of time, it can cause reduced extension speed of the device. If the device will not regain its normal function again it must be sent in for service to Össur.

WARNING: The device is not designed for use in presence of flammable gases or oxygen rich environments.

POWER SUPPLY INFORMATION

- The power supply provided with the unit is certified to comply with medical electrical equipment standards. It is only safe to use Össur recommended power supplies.
- Input: 100–240 VAC
- Input frequency: 50–60 HZ
- Input current: 0.6A
- Output voltage: 24 VDC
- Output current: 1.0 A
- Operating temperature: 0°C to 40°C (32°F to 104°F)
- Operating humidity: 5% to 95% relative humidity (non-condensing)
- Manufacturer: UE Electronic
- Model No.: UE24WCP-240100SPA

ELECTROMAGNETIC COMPATIBILITY PRECAUTIONS

SAFETY PRECAUTIONS

Because of electromagnetic compatibility (EMC) special precautions need to be taken with the device:

- The device should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the device should be checked to make sure it operates normally.
- The device can be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) phones.
- The device can be susceptible to interference from other equipment, even if the equipment complies with CISPR EMISSION requirements.

MAINTENANCE

The unit should be sent to Össur for repairs or service. No other company or individual is allowed to service or repair the unit in any way.

Before sending the device for repairs or service, make sure that the operating instructions have been followed accurately. Immediately after exposure to water the device may emit slight noise during use.

CLEANING

- Clean the knee with a damp cloth and mild soap or isopropyl alcohol when needed. Dry the product with a lint-free cloth and allow it to dry fully.
- DO NOT use solvents stronger than isopropyl alcohol.
- DO NOT use compressed air to clean the knee, as it can cause the knee to malfunction and deteriorate prematurely.

LIABILITY

The manufacturer recommends using the device only under specified conditions and for the intended purposes. The device must be maintained and used according to the instructions for use. The manufacturer is not liable for

damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the device.

CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class II device according to the classification criteria outlined in appendix IX of the guideline.

WARRANTY

For information on warranty, please see the Warranty Card

COMPLIANCE

This component has been tested according to ISO 10328 standard to three million load cycles. Depending on the amputee's activity this corresponds to a duration of use of three to five years. We recommend carrying out regular yearly safety checks

ISO 10328 - P6 - 136 kg ^{*)} 

^{*)} Body mass limit not to be exceeded!
For specific conditions and limitations of use see manufacturer's written instructions on intended use!

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby electronic equipment. The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonics Emissions IEC 61000-3-2	Class B	
Voltage Fluctuations/ Flicker Emissions vIEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, Short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms	$d=1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m	$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			where “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^o, should be less than the compliance level in each frequency range^o. Interference may occur in the vicinity of equipment marked with the following symbol: 
		Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.
		Note 2	These guidelines may not apply in situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
		*	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
		**	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1,2√P	80 MHz to 800 MHz d=1,2√P	800 MHz to 2,5 GHz d=2,3√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Type and Frequency Characteristics	Effective Radiated Power
Bluetooth 4.2 Module: Model: BT121	FID 2402 MHz~2480 MHz (1 MHz Interval, 79 Channels)	0.2663mW/MHz
BLUETOOTH 4.1 Module: Model: N5 ANT SoC	FID 2402 MHz~2480 MHz (1 MHz Interval, 79 Channels)	0.42mW/MHz

USA-FEDERAL COMMUNICATIONS COMMISSION (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.”

Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada’s website www.hc-sc.gc.ca/rpb.

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EN – Caution: Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Ossur products are used with other recommended Ossur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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