GLOBUS

DUO PRO

35 PROGRAMS

SPORT, FITNESS, BEAUTY AND HEALTH





NEW LOOK



SOFT TOUCH COVERING

Duo Pro offers a pleasant "soft touch" sensation in your hands. This special covering makes it "soft and rubbery"



VERY SIMPLE TO USE



POSSIBILITY TO PERSONALISE THE DURATION TIME OF EACH PROGRAM

High tech & Design

The new Duo Pro is a palm electro stimulator with professional features: modern and innovative design, compact dimension, back light display.

- All the best programs: in order to offer the highest performances but with a very simple use.
- Sports and performances: 8 essential programs for the training and the Sports performances





- Fitness and body: 20 programs to reach and keep a perfect physical form, the best programs for recovering the harmony and the beauty of everyone's body.
- Health and well-being: 7 specific programs to treat and eliminate the pain







35 PROGRAMS:



7 Sports Programs

Capillarisation Warm up Post competition recovery Maximal force Resistance force Explosive force Reactivity

20 Fitness Programs

Arms drainage (swollen arms) Thighs drainage Waist drainage Abdominal drainage Gluteus drainage Thighs anticellulite lipolysis Gluteus and waist anticellulite lipolysis Abdominal anticellulite lipolysis Upper limbs anticellulite lipolysis Thighs, gluteus and waist firming Abdominal firming

Upper limbs firming Upper limbs toning

Thighs, gluteus and waist toning

Abdominal toning Upper limbs sculpting

Thighs, gluteus and waist sculpting

Abdominal sculpting Face beauty Breast toning

8 Rehab Programs

Upper limbs muscular atrophy Lower limbs muscular atrophy Tens antalgic Tens endorphinic Tens anti-inflammatory (articular pains, shoulder, knee, hip) Neck pain Lumbo-sciatica Muscular contractures (tensed muscles)

ACCESSORIES

- PALM ELECTROSTIMULATOR
- BAG FOR THE TRANSPORT
- 2 CONNECTION CABLES
- 4 REUSABLE DHESIVE ELECTRODES 45X45 MM
- 4 REUSABLE DHESIVE ELECTRODES 45X90 MM
- 3 1,5 V AAA BATTERIES
- OPERATING MANUAL



TECHNICAL FEATURES

Channels

Intensity 100 lpp

Wave form Symmetrical biphasic square pulse

Frequency 1—120 Hz Amplitude $50 - 320 \mu s$

Power supply 3 1,5 V AAA batteries type LRO3

Conformity The device complies with the ISO 9001, EN 46001 standards

for electro medical devices that are MDD (Medical Device

Directive) certified.

Globus Italia srl



