



perfecta: aerosol therapy at its absolute best

Perfecta spacer mask for adults and children

- particle size 1– 5 microns
- relatively short treatment time
- greater volume for increased peripheral deposition of the drug
- easily tolerated by users of any age, thereby increasing the treatment's effectiveness

For the local (*) treatment of:

- respiratory infections of bacterial or viral origin
- asthma
- bronchitis
- pneumonia
- cystic fibrosis

(*) Localised intake of smaller quantities of drug compared to systemic treatment reduces side effects and enhances effectiveness.



perfecta
can be used with any
piston-operated aerosol device.



AIR LIQUIDE MEDICAL SYSTEMS S.p.A.
62, Via dei Prati
25073 Bovezzo (Brescia)
Tel. (+39) 030 2015911
Fax (+39) 030 2098329
www.airliquidemedicalsystems.it
www.airliquidemedicalsystems.com



As from 1st June 2009 Markos Mefar has become
Air Liquide Medical Systems S.p.A.



www.airliquide.com

Air Liquide's Healthcare business supplies medical gases & equipment to hospitals, medical services to patients at home and is active in the fight against nosocomial infections. Its 7,800 employees serve 5,000 hospitals and 300,000 patients, on five continents. Air Liquide is the world leader in gases for industry, health and the environment.



Perfecta

New spacer mask
for adults
and children



www.airliquidemedicalsystems.it
www.airliquidemedicalsystems.com



Cod. 810.223 09/2009



The **new PERFECTA spacer mask** is designed for treating disorders of the lower airways. It is specially shaped for use with adults and children. With Perfecta, **treatment is quick, safe and effective.**

Air Liquide Medical Systems has once again come up with an ingenious device offering patients of all ages a better quality of life - anywhere and at all times.



perfecta


one mask
for both adults
and children

Perfecta, which comes with the Mefar 2000 nebuliser, is a **new spacer mask** for aerosol therapy and can be used by both adults and children.

It sprays the drug **straight into the target organ**, and the whole treatment only takes **two minutes.**



Perfecta has a **larger SPACER MASK** to allow **enhanced drug deposition in the peripheral area of the lungs** and less in the mouth and pharynx.



Particle size 1-5 microns.

perfecta

is a mouth and nose mask
for children

When used with children, the mask must be **placed over the mouth and nose.**

The special shape of the mask and the soft material mean it



can also be used with **very young children**, who are not usually very cooperative.

A **valve** on the fitting relieves any overpressure.

perfecta

is a mouth and nose mask
for adults

When used by adults, the mask must **be placed over the mouth**, with the top edge up against the nose.

For effective treatment, the user should breathe in through the mouth and out through the nose.



If you breathe out through the mouth, the **valve** on the fitting relieves any overpressure.



The **PERFECTA mask** features a perfect **ergonomic design.**

It is made of silicone, **an anallergic material that is soft and pleasing to the touch.**

This makes it very easy to wear by patients of any age and improves compliance by adults and children.